Operating Instructions | Betriebsanleitung | Mode d'emploi |
Instrucciones de manejo | Manuale d'uso |
Instruções de Operação | Руководство по эксплуатации |
操作说明 | 取扱説明書 | 사용 설명서

Original Operating Instructions | Original-Betriebsanleitung |
Mode d'emploi original | Instrucciones de manejo originales |
Manuale d'uso originale | Instruções de Operação Originais |
Оригинальное руководство по эксплуатации |
原始操作说明 | 取扱説明書 (オリジナル版) | 원본 사용 설명서

YDP50

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Contents

1	About This Document	6
1.1	Validity	6
1.2	Target Groups	6
1.3	Symbols Used	6
1.3.1	Warnings in Operation Descriptions	6
1.3.2	Other Symbols Used	7
2	Safety Instructions	7
2.1	Intended Use	7
2.1.1	Modifications to the Device	8
2.2	Electrical Equipment	8
2.2.1	Damage to the Device's Electrical Equipment	8
2.2.2	Power Supply Unit and Power Plug Adapter	8
2.3	Accessories and Consumables	9
2.4	Housing Breakage	9
3	Device Description	10
3.1	Device Overview	10
3.2	Connections	11
3.3	Components for the Print Function	12
4	Installation	13
4.1	Scope of Delivery	13
4.2	Selecting an Installation Site	13
4.3	Unpacking	14
4.4	Acclimatizing	14

5	Getting Started	14
5.1	Connecting the Device to the Laboratory Device	14
5.2	Mounting or Dismounting the Power Supply Unit	15
5.2.1	Mounting the Power Supply Unit	15
5.2.2	Dismounting the Power Plug Adapter	15
5.3	Connecting the Power Supply Unit to the Power Supply	16
5.4	Switching the Device On or Off	16
5.5	Opening or Closing the Cover	17
5.6	Checking or Replacing the Ink Ribbon Cassette	17
5.6.1	Checking the Ink Ribbon Cassette	17
5.6.2	Replacing the Ink Ribbon Cassette	18
5.7	Inserting the Paper Roll into the Device	18
5.8	Establishing Communication with the Connected Laboratory Device	19
5.8.1	Connection via USB	19
5.8.2	Connection via RS232	20
6	Operation	20
6.1	Feeding Paper	20
7	Cleaning and Maintenance	21
7.1	Cleaning	21
8	Malfunctions	22
8.1	Device Malfunctions	22
8.2	Malfunctions when Printing	22
8.3	Poor Print Quality	24
9	Decommissioning	24
9.1	Decommissioning the Device	24
10	Transport	25
10.1	Transporting the Device	25
11	Storage and Shipping	25
11.1	Storing	25
11.2	Returning the Device and Parts	25

12 Disposal	26
12.1 Disposing of the Device and Parts	26
13 Technical Data	27
13.1 Dimensions and Weights	27
13.2 Ambient Conditions	27
13.3 Storage Conditions	28
13.4 Electrical Data	28
13.4.1 Power Supply.....	28
13.4.2 Safety and Electrical Compatibility.....	29
13.5 Interfaces	30
13.6 Device Functions	30
13.6.1 Printing.....	30
13.6.2 Paper.....	30
14 Accessories and Consumables	31
14.1 Accessories	31
14.2 Consumables	31
15 Sartorius Service	32
16 Conformity Documents	32
EC / EU Declaration of Conformity	295
UK Declaration of Conformity	301

1 About This Document

1.1 Validity

These instructions are part of the device; they must be read in full and stored. These instructions apply to the device in the following versions:

Device	Item Number
YDP50, printer	YDP50

1.2 Target Groups

These instructions are addressed to the following target groups. The target groups must possess the knowledge specified below.

Target Group	Knowledge and Qualifications
Operator	The operator is familiar with the device and the associated work processes. The operator understands the hazards which may arise when working with the device and knows how to prevent them.

1.3 Symbols Used

1.3.1 Warnings in Operation Descriptions

NOTICE

Denotes a hazard that may result in property damage if it is **not** avoided.

1.3.2 Other Symbols Used

- ▶ Required action: Describes activities that must be carried out. The activities in the sequence must be carried out in succession.
- ▷ Result: Describes the result of the activities carried out.
- [] Refers to operating and display elements. Text inside brackets indicates status, warning, and error messages.

2 Safety Instructions

2.1 Intended Use

The device prints data from laboratory devices, for example to print out measurement results from a laboratory balance. The device can be used in laboratories.

Printing is carried out using the “dot matrix” printing process on standard paper. The printer must be connected to the laboratory device via a cable (for suitability, see Chapter “13.5 Interfaces”, page 30). The printout can be separated from the paper roll using the paper tear bar.

The device is intended solely for use in accordance with these instructions. Any other use is considered **improper**.

Operating Conditions for the Device

Do **not** use the device in potentially explosive environments. Only use the device indoors.

The device may only be used with the equipment and under the operating conditions described in the Technical Data section of these instructions.

2.1.1 Modifications to the Device

If the device is modified: Persons may be put at risk. Device-specific documents and product approvals may lose their validity.

2.2 Electrical Equipment

2.2.1 Damage to the Device's Electrical Equipment

Damage to the device's electrical equipment, e.g., damaged insulation, may result in unforeseeable hazards.

- ▶ If the electrical equipment of the device is defective, cut off the power supply and contact Sartorius Service.
- ▶ Keep parts under voltage away from moisture. Moisture can cause short circuits.

2.2.2 Power Supply Unit and Power Plug Adapter

The use of an **unauthorized** power supply unit or power supply cable may cause life-threatening injuries as a result of electric shocks, for example.

- ▶ Only use the original power supply unit and original power plug adapter.
- ▶ If the power supply unit or power plug adapter needs to be replaced: Contact Sartorius Service. Do **not** repair or modify the power supply unit or power plug adapter.

2.3 Accessories and Consumables

The use of unsuitable accessories and consumables can affect the functionality and safety of the device and have the following consequences:

- Risk of injury to persons
 - Damage, malfunctions, or failure of the device
- ▶ Only use accessories, consumables, and spare parts that have been approved by Sartorius and are safe for use.

2.4 Housing Breakage

The plastic device may break if it falls or is handled incorrectly. Broken edges can cause cuts.

- ▶ If the device is damaged: Dispose of the device (see Chapter “12 Disposal”, page 26).

3 Device Description

3.1 Device Overview

Fig.1: YDP50 (front view)

Pos.	Name	Description
1	Cover	Covers the paper roll and the ink ribbon cassette.
2	LED display	Indicates the device operating status: <ul style="list-style-type: none"> – LED lights up: The device is ready for use. – LED flashes: There is a malfunction. – LED off: Device is switched off.
3	On Off switch	Switches the printer on or off.
4	Feed button	For moving the inserted paper forward.

3.2 Connections


Fig.2: Connections (rear view)

Pos.	Name	Description
1	RS232C	Connects the device to a laboratory device.
2	USB-B	Connects the device to a laboratory device.
3	Power supply	For connecting the supplied power supply unit.

3.3 Components for the Print Function

Fig.3: Components for the print function (top cover open)

Pos.	Name	Description
1	Roll core	Keeps the inserted paper roll in position.
2	Ink ribbon cassette	Contains the ink ribbon.
3	Paper tear bar	
4	Paper roll	

4 Installation

4.1 Scope of Delivery

Article	Quantity
Device	1
Power supply unit	1
Country-specific power plug adapter	10
Printer connection cable, type USB	1
Paper roll	2
Operating Instructions	1

4.2 Selecting an Installation Site

Procedure

- ▶ **NOTICE** The electronic equipment and the printing unit may be damaged by chemically corrosive environments, high temperatures, high humidity, and vibrations! Ensure that the installation conditions have been met (see Chapter “13.2 Ambient Conditions”, page 27).
- ▶ **NOTICE** Loss of function of the electronic equipment and the printing unit due to water spray or a dusty environment. Ensure that the device is protected against water spray and dust.

4.3 Unpacking

Procedure

- ▶ Unpack the device.
- ▶ Install the device at the intended installation site. The device must **not** be dropped or shaken.
- ▶ We recommend keeping the original packaging to return the device securely and appropriately.

4.4 Acclimatizing

When a cold device is brought into a warm environment: The temperature difference can lead to condensation from humidity in the device (moisture formation). Moisture in the device can lead to malfunctions.

Procedure

- ▶ Allow the device to acclimatize to the installation site. The device must be disconnected from the power supply during that time.

5 Getting Started

5.1 Connecting the Device to the Laboratory Device


Procedure

- ▶ **NOTICE** Improper connection may damage the device! If the device is connected to electronic components, e.g. a laboratory device: The device must be disconnected from the power supply. Ensure that the device is disconnected from the power supply.
- ▶ Connect the printer connection cable to the “USB-B” port.
- ▶ Connect the printer connection cable to the laboratory device (see laboratory device operating instructions).

5.2 Mounting or Dismounting the Power Supply Unit

5.2.1 Mounting the Power Supply Unit


Procedure

- ▶ Select the country-specific power plug adapter. The power plug adapter must be suitable for use with the wall outlet at the installation site.
- ▶ Insert the power plug adapter into the power supply unit. The grooved button must be facing forwards.
- ▶ Push the power plug adapter as far as you can until it audibly clicks into place.
- ▶ Check whether the power plug adapter is securely locked in place by pulling it gently.
- ▶ If the power plug adapter does **not** move: It is locked in place.

5.2.2 Dismounting the Power Plug Adapter

Procedure

- ▶ Press the grooved button from above and pull back on the power plug adapter.
- ▶ Push the power plug adapter out of the power supply unit and remove it.

5.3 Connecting the Power Supply Unit to the Power Supply

Procedure

- ▶ Check whether the country-specific power plug matches the power supplies at the installation site.
 - ▶ If required: Replace the country-specific power plug adapter.
- ▶ **NOTICE** Damage to the device due to incorrect input voltage! Check whether the voltage specifications on the manufacturer's ID label match those of the power supply at the installation site.
 - ▶ If the input voltage is too high or too low: Do **not** connect the device to the power supply.
 - ▶ Contact Sartorius Service.
- ▶ Connect the power supply unit to the wall outlet (supply voltage) at the installation site.

5.4 Switching the Device On or Off

Procedure

- ▶ To switch the device on: Press the On | Off switch so that it is in the "I" position.
- ▷ The LED display lights up.
- ▷ The device is ready for use.
- ▶ To switch the device off: Press the On | Off switch so that it is in the "o" position.

5.5 Opening or Closing the Cover

Procedure

- ▶ Open the cover upward and push it backward until it stops.
- ▶ To close the cover: Fold the cover downward.


5.6 Checking or Replacing the Ink Ribbon Cassette

5.6.1 Checking the Ink Ribbon Cassette

The ink ribbon cassette is installed on delivery. The ink ribbon cassette may become loose due to vibrations during transport.

Procedure

- ▶ Open the cover.


- ▶ If the ink ribbon cassette has come loose: Press the ink ribbon cassette down on both sides until it engages.
- ▶ Close the cover.

5.6.2 Replacing the Ink Ribbon Cassette


If the print is getting fuzzy or faint: Replace the ink ribbon cassette.

Procedure

- ▶ **NOTICE** Risk of damage to the device as a result of an unsuitable Ink ribbon! Check whether the ink ribbon is damaged. Only use the approved ink ribbon cassette (see Chapter “14.2 Consumables”, page 31).
- ▶ Open the cover.
- ▶ Roll back the paper roll.

- ▶ From above, push down on the right side of the ink ribbon cassette, which is marked “PUSH”.
- ▶ The left side of the ink ribbon cassette lifts up.
- ▶ Remove the ink ribbon cassette from the device.
- ▶ Insert the new ink ribbon cassette from above and press down.
- ▶ The ink ribbon cassette will engage.

- ▶ If there are creases in the ink ribbon: Turn the winding spool in the direction indicated by the arrow until the ink ribbon is taut.
- ▶ Close the cover.


5.7 Inserting the Paper Roll into the Device

Requirements

- The device is switched on.

Procedure

- ▶ **NOTICE** The device may be damaged by using non-specified paper, e.g. thermal paper! Check whether the paper is suitable for the device (for suitability, see Chapter “13.6.2 Paper”, page 30“).
- ▶ Open the cover.

- ▶ Insert the roll core (1) into the paper roll (2).
- ▶ Place the paper roll into the paper compartment.
- ▶ The end of the paper roll must point downwards.

- ▶ Insert the end of the paper roll into the slot of the printing unit.

- ▶ To feed the paper through the printing unit: Press the feed button (1) several times or press and hold until the end of the paper comes out of the printing unit.
- ▶ Close the cover.

5.8 Establishing Communication with the Connected Laboratory Device

5.8.1 Connection via USB

Communication via USB does not require any settings and works immediately after establishing the connection via USB cable.

5.8.2 Connection via RS232

If the Sartorius factory settings are set on the laboratory device to be connected, the printer will work after establishing the connection via RS232 cable.

Requirements

The laboratory device is switched on.

Procedure

- ▶ Switch off the device.
- ▶ Press and hold the feed button and switch on the device.
- ▶ Press and hold the feed button until the message "Auto-Set" is printed.
- ▶ Check the data transfer. In order to do this, trigger a print job on the laboratory device (see laboratory device instructions).
- ▷ The correct communication parameters are printed.
- ▶ If the printed content is correct: Press the feed button. This accepts the communication parameters.
- ▷ The printer returns to the normal operating mode.

6 Operation

6.1 Feeding Paper

Procedure

- ▶ To feed the paper through the printing unit: Press and hold the feed button. This pulls in the paper in a continuous process.
- ▶ To stop the process: Release the feed button.

7 Cleaning and Maintenance

7.1 Cleaning

Requirements

The printing process is complete.

NOTICE

The electronic equipment could be damaged by improper cleaning!

Liquids or dust can damage the device or the power supply unit.

- ▶ Do **not** open the power supply unit.
 - ▶ Ensure that **no** liquids or dust get into the device or the power supply unit.
-

Procedure

- ▶ Disconnect the device from the power supply.
- ▶ **NOTICE** The device could be damaged by unsuitable cleaning agents! Only use suitable cleaning agents and observe the product information for the cleaning agent used.
- ▶ Keep corrosive liquids away from the printer.
- ▶ Wipe the housing with a slightly damp cloth. In the case of more heavy contamination, use a mild soapy solution or a suitable cleaning agent.
- ▶ Only touch the power supply unit with dry hands.

8 Malfunctions

8.1 Device Malfunctions

Malfunction	Cause	Solution	Chapter, Page
The LED display is inactive.	The device is switched off.	Switch on the device.	5.4, 16
	The device is disconnected.	Connect the device to the power supply at the installation site.	5.3, 16

8.2 Malfunctions when Printing

Malfunction	Cause	Solution	Chapter, Page
The device does not print.	No paper has been inserted.	Insert a paper roll.	5.7, 18
	The ink ribbon cassette has been inserted incorrectly.	Remove the ink ribbon cassette and insert it again.	5.6.2, 18
	The device is not correctly connected to the laboratory device.	Check the connection.	5.1, 14
	The RS232 connection cable is not correctly assigned.	Use the printer connection cable, which is available as an accessory from Sartorius.	14.1, 31

Malfunction	Cause	Solution	Chapter, Page
The device does not print.	The communication parameters between the device and the laboratory device do not match.	Reset the laboratory device to the factory setting.	See the operating instructions of the laboratory device.
		Establish communication with the laboratory device.	5.8.1, 19
There is a paper jam or the printed content overlaps.	Excess paper in the device.	Open the cover and cut off the excess paper.	6.1, 20
		Press the feed button.	
	Remove the excess paper.		
	There is no paper roll core.	Insert the paper roll with the paper roll core.	5.7, 18
Noise during paper feeding	There is a foreign body in the paper feed slot.	Carefully remove the foreign body using tweezers.	

8.3 Poor Print Quality

Malfunction	Cause	Solution	Chapter, Page
The print quality is faint.	The ink ribbon is old or has worn out.	Insert a new ink ribbon cassette.	5.6.2, 18
The print quality is poor.	The ink ribbon cassette has been inserted incorrectly.	Remove the ink ribbon cassette and insert it again.	5.6.2, 18
The print quality is patchy.	The ink ribbon is wrinkled.	Remove the ink ribbon cassette.	5.6.2, 18
		Tighten the ink ribbon.	5.6.2, 18
		Insert the ink ribbon cassette.	5.6.2, 18

9 Decommissioning

9.1 Decommissioning the Device

Procedure

- ▶ Switch off the device.
- ▶ Disconnect the device from the power supply and all devices. In order to do this, remove all connection cables.
- ▶ Remove the paper roll from the device.
- ▶ Remove the ink ribbon.
- ▶ Clean the device.

10 Transport

10.1 Transporting the Device

Requirements

The device has been decommissioned.

Procedure

- ▶ Transport the device in the original packaging.

11 Storage and Shipping

11.1 Storing

Requirements

The device has been decommissioned.

Procedure

- ▶ Store the device according to the storage conditions (see Chapter "13.3 Storage Conditions", page 28).

11.2 Returning the Device and Parts

Defective devices or parts can be returned to Sartorius. Returned devices must be clean and packed in their original packaging.

Transport damage as well as measures for subsequent cleaning and disinfection of the device or parts by Sartorius are charged to the sender.

Devices contaminated with hazardous materials, e.g., harmful biological or chemical substances, will **not** be accepted for repair or disposal.

Procedure

- ▶ Decommission the device.
- ▶ Contact Sartorius Service for instructions on how to return devices or parts (see www.sartorius.com).
- ▶ Pack the device and its parts in their original packaging for return.

12 Disposal

12.1 Disposing of the Device and Parts

The device must be disposed of properly by disposal facilities.

Procedure

- ▶ Dispose of the device and packaging in accordance with local government regulations.

13 Technical Data

13.1 Dimensions and Weights

	Unit	Value
Dimensions (L x W x H)	mm	188.2 x 114.6 x 86.6
Weight incl. packaging, approx.	kg	1.8

13.2 Ambient Conditions

	Unit	Value
Installation site		
Standard laboratory rooms		
Height above sea level, maximum	m	3000
For indoor use only		
Pollution level of the power supply unit according to IEC 62368-1		2
Operating temperature	°C	0 - +40
Relative humidity		
At temperatures of up to 34°C Then linear decrease from 90% at 34°C to 65% at 40°C	%	10 - 90

13.3 Storage Conditions

	Unit	Value
Temperature	°C	- 25 - + 70
Relative humidity	%	10 - 90

13.4 Electrical Data

13.4.1 Power Supply

	Unit	Value
Power supply only permitted using the Sartorius power supply unit		
Primary		
AC voltage	V	100 - 240 (±10%)
Frequency	Hz	50 - 60 (±5%)
Secondary		
DC voltage	V	12 (±5%)
Current, maximum	A	2.0
Power supply unit fuses		
Type		Electronic
Protection class according to IEC 62368-1		II
For further data, see label on the power supply unit		

13.4.2 Safety and Electrical Compatibility

Safety according to IEC 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements

Electromagnetic compatibility, according to EN 61326-1 Electrical equipment for measurement, control, and laboratory use – EMC requirements – Part 1: General requirements (IEC 61326-1)

Suitable for use in industrial areas (interference immunity)

Interference emission, class B

Class B: Suitable for use in residential areas and areas that are directly connected to a low voltage network that supplies residential buildings.

NOTE: This equipment has been tested and found to comply with the limits for Class B digital devices, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in residential installations. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
 - Increase the distance between the equipment and receiver.
 - Connect the equipment to an outlet on a separate circuit from the one the receiver is connected to.
 - Consult your supplier or an experienced radio/TV technician for help.
-

Shielded cables: Connections between the device and peripheral equipment must be made using shielded cables in order to maintain compliance with FCC radio frequency emission limits.

Modifications: Any modifications made to this device that are not approved by Sartorius may void the authority granted to the user by the FCC to operate this equipment.

Temporary connection interruptions may occur in some cases if the printer is connected via USB-B. The printer will automatically resume all functions.

13.5 Interfaces

Type, serial
USB-B, type 2.0
RS232, D-Sub, 9-pin

13.6 Device Functions

13.6.1 Printing

	Unit	Value
Printing process: Dot matrix via 8 print magnets		
Print speed	Lines/second	1.7
Print width, maximum	mm	57
Character size (W x H)	mm	1.7 x 2.6
Supported codes: Standard ASCII (5x7)		
Service life of the printing unit	Lines	1000000

13.6.2 Paper

	Unit	Value
Type: Continuous paper		
Size of paper roll (W x ø)	mm	57 x 50

14 Accessories and Consumables

14.1 Accessories

This table contains an excerpt of accessories that can be ordered. For information on other products, please contact Sartorius.

Article	Quantity	Order Number
Printer connection cable, type RS232, 9-pin (male) on RS232, 9-pin (female)	1	YCC-D09MF
Printer connection cable, type USB-B on USB mini (OTP)	1	69Y03293
Printer connection cable, type USB-B on USB-C (in the scope of delivery)	1	YCC-USB-C-B

14.2 Consumables

This table contains a selection of the consumables that can be ordered. For information on other products, please contact Sartorius.

Article	Quantity	Order Number
Paper roll (57 x 50 mm)	10	YDP50-PAPER
Ink ribbon cassette	1	6906918

15 Sartorius Service

Sartorius Service is at your disposal for queries regarding the device. Please visit the Sartorius website (www.sartorius.com) for information about the service addresses, services provided, or to contact a local representative.

For inquiries about the system or when contacting Sartorius Service in the event of a malfunction, ensure that you have the device information, e.g., serial number, close at hand so you can pass it on to Sartorius Service. To do this, refer to the information on the manufacturer's ID label.

16 Conformity Documents

The attached documents declare the conformity of the device with the designated directives or standards.

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This table was developed according to the provisions of SJ/T 11364.

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- X: Indicates that contents of Pb, Hg, Cd, Cr⁶⁺, PBB, PBDE in at least one of the homogeneous materials used for this part is above the limit requirement in GB/T26572.
 X: 表示该有毒或有害物质至少在该部件所用的某一均质材料中的含量超出 GB/T26572 标准规定的限量要求。
 0: Indicates that contents of Pb, Hg, Cd, Cr⁶⁺, PBB, PBDE in at least one of the homogeneous materials used for this part is below the limit requirement in GB/T26572.
 0: 表示该有毒或有害物质在该部件所有均质材料中的含量均在 GB/T26572 标准规定的限量要求以下。

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SARTORIUS

EG-/EU-Konformitätserklärung
EC / EU Declaration of Conformity

Hersteller
Manufacturer Sartorius Scientific Instruments (Beijing) Co., Limited
33 Yu An Road, Shunyi District, Beijing, China

erklärt in alleiniger Verantwortung, dass das Betriebsmittel
declares under sole responsibility that the equipment

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in the form as delivered fulfils all the relevant provisions of the following European Directives and meets the applicable requirements of the harmonized European Standards including any amendments valid at the time this declaration was signed listed below:

	EMV EMC	RoHS	LVD LVD
Richtlinie Directive	2014/30/EU	2011/65/EU	2014/35/EU
Norm(en) Standard(s)	EN 61326-1:2013 EN IEC 61326-1:2021	EN IEC 63000:2018	EN IEC 62368-1:2020 +A11:2020

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Traduction du document original

SARTORIUS

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Fabricant **Sartorius Scientific Instruments (Beijing) Co., Limited**
33 Yu An Road, Shunyi District, Beijing, Chine

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Sartorius Scientific Instruments (Beijing) Co., Limited
Beijing, **16/01/2024**

Dr. Yuguang Zhao
Managing Director

Xiulin Zhu
Testing & Certificate Engineer

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SARTORIUS

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Fabricante **Sartorius Scientific Instruments (Beijing) Co., Limited**
33 Yu An Road, Shunyi District, Pekín, China

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SARTORIUS

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Fabbricante **Sartorius Scientific Instruments (Beijing) Co., Limited**
33 Yu An Road, Shunyi District, Beijing, Cina

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SARTORIUS

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Fabricante **Sartorius Scientific Instruments (Beijing) Co., Limited**
33 Yu An Road, Shunyi District, Beijing, China

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Original

UK Declaration of Conformity

Manufacturer Sartorius Scientific Instruments (Beijing) Co., Limited
33 Yu An Road, Shunyi District, Beijing, China
declares under sole responsibility that the equipment

Device type Printer

Model YDP50

in the form as delivered fulfils all the relevant provisions of the following British Regulations and meets the applicable requirements of the British Designated Standards including any amendments valid at the time this declaration was signed listed below:

	EMC	RoHS	LVD
Regulation	The Electromagnetic Compatibility Regulations 2016 UK Statutory Instruments 2016 No. 1091	The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012 UK Statutory Instruments 2012 No. 3032	The Electrical Equipment (Safety) Regulations 2016 UK Statutory Instruments 2016 No. 1101
Standard(s)	BS EN 61326-1:2013 BS EN IEC 61326-1:2021	BS EN IEC 63000:2018	BS EN IEC 62368-1 +A11:2020

The person authorised to compile the technical file:

Sartorius UK Ltd.
Longmead Business Centre, Blenheim Road
KT19 9 QQ Epsom, Surrey, UK

Sartorius Scientific Instruments (Beijing) Co., Limited
Beijing, 2024-01-16

Dr. Yuguang Zhao
Managing Director

Xiulin Zhu
Testing & certificate engineer

Sartorius Scientific Instruments (Beijing)
Co., Limited
No.33 Yu'an road, Airport Economic Zone,
Shunyi District, 101300 Beijing, P.R. China

Phone: +86 10 8042 6300 +86 800 8100 124
www.sartorius.com

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote all genders.

Copyright notice:

These instructions, including all components, are protected by copyright.

Any use beyond the limits of the copyright law is not permitted without our approval.

This applies in particular to reprinting, translation and editing irrespective of the type of media used.

Last updated:

01 | 2024

© 2024
Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Str. 20
37079 Goettingen, Germany

JW | Publication No.: WYD6165ep240102