
Operating Instructions | Betriebsanleitung | Mode d'emploi | Manuale di istruzioni |
Instrucciones de manejo | Инструкция по эксплуатации | 操作手册 | 取扱説明書

YDP30

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The image shows a printer's self-test printout with various fields and test patterns. Lines connect specific fields to their descriptions:

- PRINTER INFO. (Printer model and main board firmware version)
- YDP30 Version: 6.92 E2 (Printer model name and main board firmware version)
- SERIAL NO. (Serial number)
- MILAGE(m): 661 (Print distance)
- CHECKSUM: 0783361C TCF (Main board firmware checksum)
- SERIAL PORT: 9600,0,8,1 (Serial port settings)
- CODE PAGE: 850 (Code page)
- COUNTRY CODE: 001 (Country code)
- SPEED: 2 INCH (Print speed)
- DENSITY: 8.0 (Print density)
- SIZE: 1.66 1.58 (Label size (width, height))
- GAP: 0.00 0.00 (Gap size (vertical gap, offset))
- TRANSPARENCY: 7 (Sensor sensitivity)
- ***** (Separator line)
- FILE LIST: (File management information)
- DRAM FILE: 0 FILE(S)
- FLASH FILE: 2 FILE(S)
- LOGO.BMP 22246 BYTES
- AUTO.BAS 2951 BYTES
- PHYSICAL DRAM: 8192 KBYTES
- AVAILABLE DRAM: 256 KBYTES FREE
- PHYSICAL FLASH: 4096 KBYTES
- AVAILABLE FLASH: 2535 KBYTES FREE
- END OF FILE LIST
- ***** (Separator line)
- [Test patterns]

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in the form as delivered complies with the essential requirements of the following European Directives and meets the applicable requirements of the harmonized European Standards listed below:

2004/108/EG
2004/108/EC Elektromagnetische Verträglichkeit
Electromagnetic compatibility
EN 61326-1:2006
Elektrische Mess-, Steuer-, Regel- und Laborgeräte - EMV- Anforderungen – Teil 1: Allgemeine Anforderungen
Electrical equipment for measurement, control and laboratory use - EMC requirements – Part 1: General requirements

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Sartorius Lab Instruments GmbH & Co. KG
Goettingen, 2013-12-19

Dr. Reinhard Baumfalk
Vice President R&D

Dr. Dieter Klausgrete
Head of International Certification Management

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This declaration certifies conformity with the above mentioned EC and EU Directives, but does not guarantee product attributes. Unauthorised product modifications make this declaration invalid. The safety information in the associated product documentation must be observed.

Sartorius Lab Instruments GmbH & Co. KG
Weender Landstrasse 94-108
37075 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius.com

Copyright by
Sartorius Lab Instruments GmbH & Co. KG,
Goettingen, Federal Republic of Germany
No part of this publication may be
reprinted or translated in any form or by
any means without prior written
permission from Sartorius.
All rights reserved by Sartorius in
accordance with copyright law.
The information and figures contained
in these instructions correspond to the
version date specified below.
Sartorius reserves the right to make
changes to the technology, features,
specifications, and design of the
equipment without notice.

Date:
April 2014,
Sartorius Lab Instruments GmbH & Co. KG,
Goettingen, Germany