Instructions for Use | Gebrauchsanleitung  
Manuel d'utilisation | Instrucciones de uso | Istruzioni per l'uso  
使用说明 | 使用説明書 | 사용 지침 | שומיש תוארה

## Vivaflow<sup>®</sup> SU

Single Use TFF Cassettes | Single-use-TFF-Kassetten  
Cassettes TFF à usage unique | Cartuchos TFF de un solo uso  
Cassette monouso per TFF | 一次性 TFF 膜包 | シングルユースTFFカセット  
일회용 TFF 카세트 | מעפ דח שומישל TFF תוינסחמ

ZUNT-No.: 1000130516

SARTORIUS

English	page	3
Deutsch	Seite	24
Français	page	47
Español	página	70
Italiano	pagina	93
中文	页码	116
日本語	ページ	136
한국어	페이지	157
עברית	עמוד	178

# Content

<b>1</b>	<b>About these Instructions</b>	<b>5</b>
1.1	Validity	5
1.2	Related Documents	6
1.3	Target Groups	6
1.4	Symbols Used	6
1.4.1	Warnings in Operation Descriptions	6
1.4.2	Other Symbols Used	6
<b>2</b>	<b>Safety Instructions</b>	<b>7</b>
2.1	Intended Use	7
2.1.1	Operating Conditions for the Product	7
2.1.2	Modifications to the Product	7
2.2	Qualifications of Personnel	7
2.3	Functionality of the Product Parts	8
2.4	Personal Protective Equipment	8
2.5	Accessories	8
2.6	Chemical Compatibility	8
<b>3</b>	<b>Product Description</b>	<b>9</b>
3.1	Product Overview	9
3.2	Filtration Process	10
<b>4</b>	<b>Process Preparation</b>	<b>11</b>
4.1	Scope of Delivery	11
4.2	Unpacking	11
4.3	Connecting Vivaflow® Cassettes (optional)	12
4.4	Setting Up the Cassettes for Ultrafiltration   Dialfiltration	13
4.5	Pre-Washing the Cassette(s) (optional)	16

<b>5</b>	<b>Application</b> .....	<b>17</b>
5.1	Performing Ultrafiltration.....	17
5.2	Performing Diafiltration.....	18
5.3	Collecting Retentate .....	19
5.4	Malfunctions.....	19
<b>6</b>	<b>Technical Data</b> .....	<b>21</b>
6.1	Dimensions and Weights .....	21
6.2	Materials Vivaflow® Cassette and Pressure Indicator .....	21
6.3	Ambient Conditions .....	22
6.4	Operating Conditions .....	22
6.5	Shelf Life.....	23
<b>7</b>	<b>Accessories</b> .....	<b>23</b>

# 1 About these Instructions

## 1.1 Validity

These instructions are part of the product. These instructions apply to the following versions of the product:

Vivaflow® SU	Type	
	2pc, individual operation	2pc, serial operation
PES MWCO		
5 kDa	VF-S050P0005-IV	VF-S050P0005-SV
10 kDa	VF-S050P0010-IV	VF-S050P0010-SV
30 kDa	VF-S050P0030-IV	VF-S050P0030-SV
50 kDa	VF-S050P0050-IV	VF-S050P0050-SV
100 kDa	VF-S050P0100-IV	VF-S050P0100-SV
300 kDa	VF-S050P0300-IV	VF-S050P0300-SV
1,000 kDa	VF-S050P1000-IV	VF-S050P1000-SV
0.2 µm	VF-S050P2000-IV	VF-S050P2000-SV
RC MWCO		
2 kDa	VF-S050H0002-IV	VF-S050H0002-SV
5 kDa	VF-S050H0005-IV	VF-S050H0005-SV
10 kDa	VF-S050H0010-IV	VF-S050H0010-SV
30 kDa	VF-S050H0030-IV	VF-S050H0030-SV
100 kDa	VF-S050H0100-IV	VF-S050H0100-SV
300 kDa	VF-S050H0300-IV	VF-S050H0300-SV

## 1.2 Related Documents

- ▶ In addition to these instructions, please observe the following documents:
  - Instructions for the peristaltic pump drive
  - Instructions for the peristaltic pump head

## 1.3 Target Groups

These instructions are addressed to the following target groups. The target groups must possess the knowledge specified below.

Target group	Knowledge and Qualifications
Operator	The operator is familiar with the product and the associated work processes. The operator understands the hazards which may arise when working with the product, and knows how to prevent them.

## 1.4 Symbols Used

### 1.4.1 Warnings in Operation Descriptions

#### **WARNING**

Denotes a hazard that may result in death or severe injury if it is **not** avoided.

### 1.4.2 Other Symbols Used

- ▶ Required action: Describes activities that must be carried out. The activities in the sequence must be carried out in succession.
- ▷ Result: Describes the result of the activities carried out.

## Figures in These Instructions

Depending on the product configuration, the figures depicting the product may differ slightly from the supplied product. The variants shown in these instructions are examples.

## 2 Safety Instructions

### 2.1 Intended Use

Vivaflow® SU is a ready-to-use tangential flow filtration cassette designed for ultrafiltration and | or diafiltration of macromolecules in biological and aqueous solutions.

The product is supplied non-sterile. It is intended for single use and must be disposed of after one use.

The product may only be used for general laboratory use.

The product is intended exclusively for use in accordance with these instructions. Any other use is considered **improper**.

#### 2.1.1 Operating Conditions for the Product

The product may only be used with the equipment and under the operating conditions described in the Technical Data section of these instructions.

#### 2.1.2 Modifications to the Product

If the product is modified: Persons may be put at risk. Product-specific documents and product approvals may lose their validity.

For queries regarding modifications to the product, contact Sartorius.

### 2.2 Qualifications of Personnel

Persons who do **not** possess adequate knowledge about how to use the product safely may injure themselves and other persons.

If a particular qualification is required for an activity: The target group will be specified. If **no** qualification is specified: The activity may be carried out by the “operator” target group.

## 2.3 Functionality of the Product Parts

**Non-functioning** product parts, e.g., as a result of damage or wear, can cause malfunctions. There is a risk of injury to persons.

- ▶ If product parts are **not** functional: Do **not** use the product.

## 2.4 Personal Protective Equipment

Personal protective equipment protects against risks arising from the product. If the personal protective equipment is missing or is unsuitable for the work processes on the product: Persons may be injured.

The following personal protective equipment must be worn:

- Protective work clothing
- Safety gloves
- Safety glasses

## 2.5 Accessories

The use of unsuitable accessories can affect the functionality and safety of the product and have the following consequences:

- Risk of injury to persons
  - Damage, malfunctions, or failure of the product
- ▶ Only use accessories that have been approved by Sartorius and are safe for use.

## 2.6 Chemical Compatibility

The solutions used in the filtration process must be compatible with the materials of the product. Unsuitable solutions can cause damage to the product, e.g. leakage. This may result in people coming into contact with the solutions and people can be injured.

- ▶ Check the compatibility of the solutions used with the product materials (see chapter 6.2, page 21).

# 3 Product Description

## 3.1 Product Overview


Fig.1: Vivaflow® SU, front side

Pos.	Name
1	Feed port
2	Product label (not visible)
3	Retentate port
4	Interlock (for assembly onto stand and/or serial operation)
5	Permeate port

## 3.2 Filtration Process

The filtration process can be performed when the cassette is connected to a peristaltic pump. Macromolecules sufficiently larger than the nominal pore size of the membrane are retained in the flow path and recirculated between a feed reservoir and the cassette, while solvent and microsolute permeate the membrane. This process results in progressive concentration of the macromolecule of interest. With the addition of a suitable exchange buffer, the product can also be used in diafiltration mode for sample desalting or buffer exchange.

# 4 Process Preparation

## 4.1 Scope of Delivery

Article	Quantity
Vivaflow® SU cassette	2
Tubing kit (for individual or serial operation)	1
Quick start guide	1

## 4.2 Unpacking

### Procedure

- ▶ Check that the packaging material is not damaged.
- ▶ Check the expiry date on the packaging label.
- ▶ Open the packaging.
- ▶ Remove and discard the caps from the feed, permeate and retentate ports.


### 4.3 Connecting Vivaflow® Cassettes (optional)

For processing volumes up to 1L or for faster ultrafiltration | diafiltration, a maximum of 2 cassettes can be operated in series.

Material: Cassette stand

#### Procedure

- ▶ Slide the interlock on the base of one cassette into the interlock on the top of another cassette.
- ▶ Connect the retentate port from the lower cassette to the feed port on the upper cassette, using the serial connector.
- ▶ Slide the interlock on the base of the lower cassette into the cassette stand.

## 4.4 Setting Up the Cassettes for Ultrafiltration | Dialfiltration


Fig. 2: Setup for ultrafiltration

Pos.	Name
1	Retentate tube
2	Feed tube
3	Feed reservoir
4	Peristaltic pump head
5	Peristaltic pump drive
6	Pressure indicator
7	Vivaflow® cassette
8	Permeate vessel
9	Permeate tube
10	Cassette stand


Fig.3: Setup for diafiltration

Pos.	Name
1	Feed reservoir
2	Retentate tube
3	Feed tube
4	Exchange buffer tube
5	Exchange buffer vessel
6	Peristaltic pump head
7	Peristaltic pump drive
8	Pressure indicator
9	Vivaflow® cassette

Pos.	Name
10	Permeate vessel
11	Permeate tube
12	Cassette stand

Material:   – Cassette stand  
               – Feed and permeate vessels

### Requirements

The peristaltic pump head is installed on the peristaltic pump drive (see instructions for the peristaltic pump head and drive).

### Procedure

- ▶ Slide the interlock on the base of the cassette into the cassette stand.
- ▶ Connect the feed tube to a suitable feed vessel, e.g. feed reservoir, and to the feed port of the cassette.
  - ▶ If two cassettes will be operated in series: Connect the feed tube to the feed port on the lower cassette.
- ▶ Insert the feed tube into the peristaltic pump head (see instructions for the peristaltic pump head).
- ▶ Connect the permeate tube(s) to the permeate port(s) and place it | them into a suitable permeate vessel, e.g. laboratory bottle.
- ▶ Connect the pressure indicator to the retentate port of the cassette.
  - ▶ If two cassettes will be operated in series: Connect the pressure indicator to the retentate port on the upper cassette.
- ▶ Connect the retentate tube to the pressure indicator on the cassette and to the feed reservoir.

## 4.5 Pre-Washing the Cassette(s) (optional)

Pre-Washing the cassette(s) can be useful to confirm that the tubing connections are secure before starting ultrafiltration | diafiltration.

Material: 100-200 mL water or buffer

### Procedure

- ▶ Place 100 - 200 mL deionized water or buffer in the feed reservoir.
- ▶ Start the peristaltic pump with a suitable flow rate (see chapter “6.4 Operating Conditions”, page 22).
- ▷ The flushing process starts.
- ▷ The pressure indicator should read up to 2.5 bar.
- ▶ Check the system for leaks.
- ▶ When approximately 50 - 100 mL of deionized water or buffer is in the permeate vessel: Stop the peristaltic pump.
- ▷ The flushing process is complete.
- ▶ Empty the feed reservoir and the permeate vessel.

# 5 Application

## 5.1 Performing Ultrafiltration

### Requirements

- The cassette(s) is | are set up for ultrafiltration.
- The cassette(s) has | have been flushed (optional).

---

### WARNING

#### **Risk of infection due to leakage resulting in contact with solutions!**

- ▶ During the process, check the system for leaks.
  - ▶ Observe the operating conditions (chapter “6.4 Operating Conditions”, page 22).
- 

### Procedure

- ▶ Place the solution in the feed reservoir.
- ▶ Start the peristaltic pump with a suitable flow rate (see chapter “6.4 Operating Conditions”, page 22).
- ▷ The ultrafiltration process starts.
- ▷ The pressure indicator should read up to 2.5 bar.
- ▶ When the desired retentate or permeate volume is reached:
  - ▶ Reduce the flow rate to 20-40 mL/min.
  - ▶ Circulate the sample for another 1-2 minutes.
- ▶ Stop the peristaltic pump.
- ▷ The ultrafiltration process is complete.

## 5.2 Performing Diafiltration

- Material:
- 5x sample volumes of exchange buffer
  - Diafiltration tubing
  - Exchange buffer vessel

### Requirements

The sample volume is  $\leq 500$  mL.

---

### WARNING

#### Risk of infection due to leakage resulting in contact with solutions!

- ▶ During the process, check the system for leaks.
  - ▶ Observe the operating conditions (chapter “6.4 Operating Conditions”, page 22).
- 

### Procedure

- ▶ Empty the permeate vessel.
- ▶ Fill the exchange buffer vessel, e.g. laboratory bottle, with 5x sample volumes of an appropriate exchange buffer.
- ▶ Connect exchange buffer tubing to the exchange port on the feed reservoir and place the other end to the bottom of the exchange buffer vessel.
- ▶ Start the peristaltic pump with a suitable flow rate (see chapter “6.4 Operating Conditions”, page 22).
- ▷ The diafiltration process starts.
- ▷ The pressure indicator should read up to 2.5 bar.
- ▶ When the required volume of exchange buffer has been added to the system and the desired retentate volume has been reached: Stop the peristaltic pump.
- ▷ The diafiltration process is complete.

## 5.3 Collecting Retentate

Material: 25-50 mL water or buffer

### Procedure

- ▶ Disconnect the feed tube from the feed reservoir.
- ▶ Start the peristaltic pump.
- ▷ The residual retentate in the cassette(s) and tubing is pumped back into the feed reservoir.
- ▶ Stop the peristaltic pump.
- ▶ For a more thorough recovery, flush the system:
  - ▶ Fill 25-50 mL of water or sample buffer into a suitable vessel.
  - ▶ Place the feed tube into the vessel.
  - ▶ Start the peristaltic pump.
  - ▷ The buffer flows through the cassette(s) and tubing and is collected into the retentate, in the feed reservoir.
  - ▶ Stop the peristaltic pump.

## 5.4 Malfunctions

Fault	Cause	Solution
Liquids leak out of the system during flushing   filtration	The tubes are <b>not</b> connected properly.	Check and tighten the tubing connections.
Filtration does <b>not</b> occur.	The wrong tube is inserted in the peristaltic pump head.	Ensure that the feed tube is inserted into the peristaltic pump head.
	The feed tube is over- or under-compressed.	Use a peristaltic pump head that is compatible with the feed tube dimensions.

Fault	Cause	Solution
Filtration does <b>not</b> occur.	The flow rate is too low.	If the system pressure is lower than the operating limit, increase the flow rate on the peristaltic pump. <hr/>
		If the flow rate cannot be increased, check that the peristaltic pump and pump head can deliver the recommended flow rate. <hr/>

## 6 Technical Data

### 6.1 Dimensions and Weights

	Unit	Value
Cassette dimensions, excluding ports (L x W x H)	mm	98 x 13 x 116
Tubing dimensions (L x ID x WT)		
Feed tube	mm	1,000 x 3.2 x 1.6
Retentate tube	mm	1,000 x 3.2 x 1.6
Permeate tube	mm	1,000 x 1.8 x 0.6
Cassette weight	g	130

### 6.2 Materials Vivaflow® Cassette and Pressure Indicator

Component	Material
Fittings	Polyamide (PA) Polypropylene (PP)
Gasket	Silicone (SIL)
Housing	Polycarbonate (PC)
Membrane support	Polyethylene (HDPE)
Pressure indicator	Polyamide (PA) Polyoxymethylene (POM) Polypropylene (PP) Silicone (SIL) Stainless steel (SS)
Tubing	Polyvinyl chloride (PVC)

### 6.3 Ambient Conditions

	Unit	Value
Storage		
Temperature	°C	+4 - +30
Keep dry		
Keep away from light		
Operation		
Temperature	°C	+4 - +40

### 6.4 Operating Conditions

	Unit	Value
Throughput per cassette, recommended		
Minimum	mL	100
Maximum	mL	500
Flow rate, recommended	mL/min	200 - 400
Pressure, maximum		
At feed port	bar	3
At retentate port	bar	2.5
Feed tube run time, maximum*	h	6
Operating temperature	C°	+4 - +40
pH compatibility	pH	4 - 9

\* For longer run times, re-position the feed tubing in the pump head.

## 6.5 Shelf Life

	Unit	Value
Shelf life	years	2

# 7 Accessories

This table contains an excerpt of the accessories that can be ordered. For information on other products, contact Sartorius.

Item	Quantity	Order number
Peristaltic Pump Drive	1	VF-APD0001-1
Peristaltic Pump Head	1	VF-APH0001-1
Cassette Stand	1	VFA016
Diafiltration tubing	1	VF-ATD0001-1
Feed Reservoir	1	VFA006

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Sartorius Stedim Lab Ltd.  
Sperry Way, Stonehouse  
GL10 3UT, UK

Phone: +44 1453 821972  
[www.sartorius.com](http://www.sartorius.com)

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote all genders.

Copyright notice:

These instructions, including all components, are protected by copyright.

Any use beyond the limits of the copyright law is not permitted without our approval.

This applies in particular to reprinting, translation and editing irrespective of the type of media used.

Last updated:

06 | 2024

© 2024

Sartorius Stedim Lab Ltd.  
Sperry Way, Stonehouse  
GL10 3UT, UK

ML | Publication No.: SLU6135ha240601